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

## UroLift® System Deployment Sequence

- 1 Safety lock released
- 2 Needle trigger fired
  - Deploys needle and capsular tab to the capsular side of the prostate
- 3 Needle retraction lever retracted
  - Delivers capsular tab with suture under tension
- 4 Urethral release pressed
  - Attaches urethral end-piece and cuts excess suture

## Implant Technique

- Use the distal sphincter as a fulcrum
- Create anterior channel
  - Bladder neck visible from the veru
- Delivery handle parallel to the floor
- Device tip in anterior 1/3 of the urethra for lateral lobe implants

5

**Teleflex**  
INTERVENTIONAL UROLOGY

## Implant Sequence

1. Treat proximal prostate  
( $>1.5$  cm distal to bladder neck)
2. Treat distal prostate  
(anterior to verumontanum)
3. If channel not continuous, treat prostate  
between openings or as needed
4. Treat median lobe if it remains obstructive

Confirm urethral end-piece is opposed to mucosa and implants are not extending into the bladder

6

**Teleflex**  
INTERVENTIONAL UROLOGY

## Extensive Clinical Evidence Supports the UroLift System

<b>2 RCTs</b>	<b>12,000+ patients studied</b>	<b>145+ peer-reviewed clinical publications*</b>	<b>Long-term data in controlled and real-world settings</b>	<b>Healthcare claims and utilization analysis (presented at AUA 2021)</b>
---------------	---------------------------------	--	---	---

\*Select sponsored publications represented on the timeline

## The UroLift Implant Induces Tissue Remodeling

### Long-Term Biological Mechanism of Action

#### Histology of canine prostate section with implant

Compression leads to ischemia (reduced blood flow)

Moderate lobular atrophy and mild chronic inflammation in tissue surrounding implant

Ischemia causes tissue atrophy (shrinkage of tissue)

Tissue remodeling is induced by localized compression between the capsular tab and urethral end piece.

Encapsulation of implant may occur as early as 6 months.

\*The UroLift implant is customized in situ for unique anatomical variations of the prostate. Roehrborn et al., Prostate Can Prost Dis. 2021. Study sponsored by Teleflex Incorporated or its affiliates.

## “Large clinical need for an effective treatment ... less invasive than surgery” – AUA<sup>1</sup>

### From healthy bladder to permanent damage

Healthy bladder

Bladder worsens

Permanently damaged

“Since many men discontinue medical therapy, yet proportionately few seek surgery, there is a large clinical need for an effective treatment that is less invasive than surgery. With this treatment class, perhaps a significant portion of men with BOO who have stopped medical therapy can be treated prior to impending bladder dysfunction.”

American Urological Association


1. AUA Guideline 2021

## Opportunity to Change BPH Care to Better Reflect What Patients Really Want

Men's Perceptions of Importance When Considering a Minimally Invasive Procedure for BPH

10 = Most Important

### A Balancing Act

Patient priorities in BPH care

Physician discretion in BPH workup and treatment

Source: Teleflex-sponsored Men's Health Survey, November 2016, n = 1,000 men (66% of respondents were over age 45 years of age). Respondents compared items to each other by placing them in order of importance. IIC00701 -01 Rev 0.

## UroLift® System: Real Durability, Real Data Integrity, Real Patient Experience

BPH Treatment	Durability	Patient Experience			
	Surgical Retreatment Rate	Catheterization	Grade 3+ or Serious Adverse Events	Sexual Function	Hospital Stay
UroLift® System	13.6% at 5 years <sup>1</sup> Loss to Follow-up rate <sup>1</sup> : 12.8%	Lowest catheter rate of leading BPH procedures (0.7 – 0.9 days mean duration <sup>1</sup> )	0.7% <sup>1</sup>	0% EJD <sup>1,2</sup>	0 days <sup>1</sup>
Rezūm™	4.4% at 5 years <sup>2</sup> Total/Loss to Follow-up rate <sup>2</sup> : 22% (13.2% Loss to Follow-up; 8.8% decided not to participate in long-term follow-up)	>75% catheterized; 3.4 days mean duration <sup>2</sup>	1.5% <sup>2</sup>	3.0% EJD <sup>2,3</sup>	0 days <sup>2</sup>
TURP	7.6% at 2 years <sup>3</sup>	100% catheter for 1-day median <sup>3</sup>	7.6% <sup>3,4</sup>	36% EJD <sup>3,5</sup>	1.4 days <sup>3,4</sup>
Aquablation <sup>6</sup>	5.2% at 5 years <sup>4</sup> Total/Loss to Follow-up rate <sup>4</sup> : 42.1% (8.6% Loss to Follow-up; 34.2% not enrolled in long-term follow-up)	100% catheter for 1-day median <sup>4,5</sup>	6.9% <sup>4,5</sup>	10% EJD <sup>4,5</sup>	1.4 days <sup>4,5</sup>

1. Roehrborn, Can J Urol 2017; 2. McVary, J Urol 2021; 3. Bachmann, Eur Urol 2016; 4. Gillig, Can J Urol 2022; 5. Gillig, J Urol 2016; 6. McVary, J Urol 2015

<sup>1</sup>Per the FDA, Loss to Follow Up (LTFU) is defined as the act of concluding participation, prior to completion of all protocol-required elements, in a trial by an enrolled subject. The Total LTFU rate includes subjects disclosed as LTFU as well as those not enrolled or who decided not to participate in long-term follow-up. <https://www.fda.gov/oc/100275>

<sup>2</sup>De novo ED or EJD among sexually active men

<sup>3</sup>Baseline erectile and ejaculatory function was not assessed. MAC022101-01 Rev A

**Teleflex**  
INTERVENTIONAL UROLOGY

# UROLIFT

## CELEBRATING

# 450,000

## MEN TREATED FOR BPH WITH THE UROLIFT™ SYSTEM\*

An estimated 2.5 million UroLift™ implants have been delivered to patients.\*

\*Management estimates based on product sales as of September 2023. Data on file Teleflex Interventional Urology.  
Rx Only. Visit [www.teleflexurology.com](http://www.teleflexurology.com)

©2023 Teleflex Incorporated. All rights reserved. MAC022101-01 Rev F

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- Blood and/or urine – sign TMU-BB inform consent, collect blood or urine, send consent and samples to TMUH-BB RA
- Tumor and/or adjacent normal control – sign TMU-BB inform consent at OPD clinic or admission ward, send consent to operation room, call RA to OR, send consent and fresh surgical specimen to department of pathology by RA, collect tissue by **pathologist**, send tissue to TMU-BB by RA

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4

衛生福利部雙和醫院  
(委託臺北醫學大學興建經營)  
Taipei Medical University - Shuang Ho Hospital,  
Ministry of Health and Welfare

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The image shows two forms from the National Tsing Hua University Hospital. The left form is a '同意書' (Consent Form) for specimen collection, and the right form is an '申請單' (Application Form). Both forms contain fields for patient information, specimen details, and signatures.

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The image shows a '標本送出單' (Specimen Transfer Form) with a grid layout. It contains handwritten entries in blue ink, including patient names, specimen types, and dates. The form is used to track the collection and transfer of specimens.

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The image shows a '檢體簽收單' (Specimen Receipt Form) with a grid layout. It contains handwritten entries in blue ink, including patient names, specimen types, and dates. The form is used to track the receipt and processing of specimens.